

MultiFlo™ FX Multi-Mode Dispenser

MultiFlo™ FX Multi-Mode Dispenser is a versatile tool for liquid handling workflows. MultiFlo FX automates fast dispensing and washing, gentle media exchange for non- or loosely-adherent cell based assays, and dispensing into individual wells. Its unique Parallel Dispense™ design allows up to four independent reagents to be dispensed in parallel without cross-contamination.


MultiFlo FX can incorporate one or two peristaltic dispense pumps, two syringe dispense pumps, plus an available microplate wash module. The unique RAD™ (Random Access

Dispense) module can automatically dispense varying volumes into discrete wells of a 96- or 384-well plate for normalization protocols. The new patent-pending AMX™ (Automated Media Exchange) module gently automates critical steps in spheroid and non-adherent cell based assays – all one on compact platform.


MultiFlo FX, when integrated to the BioSpa 8 Automated Incubator with a BioTek imager or multi-mode reader provides complete walkaway workflow automation for many biochemical and cellular assays.



Wash 6- to 384-well plates with the wash module.



Dispense to custom-mapped wells with RAD.



AMX module enables gentle media exchange for 3D cell cultures.



Features:

- Multi-mode dispensing replaces up to four dispensers and a washer
- AMX™ Automated Media Exchange module protects 3D cell structures (e.g. spheroids, tumoroids), and suspension cells
- RAD™ Random Access Dispense enables mapped dispensing to individual wells for normalization applications
- Wash module for 6- to 384-well plates combines dispensing and washing automated on one platform
- Parallel Dispense technology: peristaltic or syringe dispensing of up to four reagents with no cross-contamination
- Cell-friendly angled dispense and wash tubes, adjustable flow rates using lidded plates
- BioSpa™ 8 Automated Incubator compatible for live and fixed cell assay automation
- Compatible with Agilent BenchCel™ Microplate Handler to automate many workflows

Typical Applications:

- Primary/secondary screening assays
- Compound storage
- Genomics and proteomics research
- ELISA
- Cell-based washing, fixing and staining
- Volume or concentration normalization


Configurations:

- MFXP1: MultiFlo FX with 1 module arm, 1 peristaltic dispense pump
MFXP2: MultiFlo FX with 2 module arms, 1 peristaltic dispense pump
MFXP2R: MultiFlo FX with 2 module arms and RAD module

See website or price list for complete listing.

Optional Modules:

- Wash module
- AMX module
- RAD module
- Dual syringe pump module
- Secondary peri-pump module



AMX enables gentle media exchange for 3D cell structures.

Optional Accessories:

- BioStack 4 Microplate Stacker
- BioSpa™ 8 Automated Incubator
- Liquid Handling Control™ Software
- Product Qualification Package
- Agilent BenchCel™ Microplate Handler



MultiFlo™ FX integrates with BioStack™ Microplate Stacker for automated processing of up to 50 plates.

Technical Details:

General

Microplate types:	Dispense, wash, RAD: 6- to 384-well, AMX: 96- and 384-well (manifold dependent)
User Interface:	Color touchscreen. Create, edit or run multiple protocols onboard.
Software:	LHC2 Software LHC2 Secure for 21 CFR Part 11 compliance (option) SiLA Compliant driver (option)
Shaking and soaking:	Programmable up to 60 minutes
Automation:	BioStack and 3rd party automation compatible BioSpa 8 Automated Incubator compatible

Dispensing: Peristaltic Pump (Multi-Channel)

Fluid delivery:	1 or 2 peristaltic pumps
Dispense speed:	96 wells, 5 µL cassette, 10 µL/ well: 3 seconds 384 wells, 1 µL cassette, 1 µL/ well: 6 seconds 1536 wells, 1 µL cassette, 1 µL/ well: 21 seconds
Dispense volume range:	500 nL - 3,000 µL/well, selectable in 1 µL increments
Flow rates:	User programmable rates from high to low
Dispense performance:	<u>1 µL cassette:</u> recommended range: 1 - 50 µL Accuracy: +5% at 1 µL, Precision: <5% CV at 1 µL <u>5 µL cassette:</u> recommended range: 5 - 2,500 µL Accuracy: +2.0% at 5 µL, Precision: <2.5% CV at 5 µL <u>10 µL cassette:</u> recommended range: 10 - 3,000 µL Accuracy: +2.0% at 10 µL, Precision: <2.0% CV at 10 µL

Dispensing: Syringe Pump (Multi-Channel)

Dispense speed:	20 µL /well, 1 x16 tubes, 96/384: 5 s/14 s 3 µL /well, 1536 wells, 2 x 32 tubes: 7 seconds
Volume range:	3 - 3,000 µL/well selectable in 1 µL increments
Dispense accuracy:	±1 µL at 5 µL and 20 µL; ±1 % at 100 µL
Dispense precision:	<2.5% CV at 20 µL; <1% CV at 100 µL

Washing

Wash volume range:	20 - 30,000 µL/well
Wash cycles:	1-10
Wash speed:	96 wells, 8-tube manifold, 3 cycles, 300 µL/well: <130 seconds
Dispense accuracy:	±3%
Dispense precision:	96-/384-well plates, 300 µL/well: <3% CV 6-well plates, 5560 µL/well: <5% CV
Residual volume:	96-well plate, 300 µL/well: <2 µL/well
Flow rates:	140 - 422 µL/well
Supply/waste bottles:	2 L, waste bottle level detection

Media Exchange: AMX (Automated Media Exchange Module)

Manifold types:	Two 8-channel autoclavable manifolds
Cassettes:	Autoclavable cassettes with 5 µL tubing
Performance:	Precision: ≤5% CV, Accuracy: ≤5%
Aspiration uniformity:	≤5%

Dispensing: RAD (Random Access Module)

Other labware:	96-well cluster tubes (minitubes) up to 50 mm height (requires custom carrier)
Manifold types:	RAD single, with plastic or steel tip with 1, 5 or 10 µL tubing, 7° angle RAD 8-to-1 plastic tip, with 5 µL tubing, angled bulk dispense chute
Volume range:	500 nL - 30,000 µL
Minimum prime volume:	1 µL cass, 18": 90 µL ; 1 µL cass, 30": 150 µL 5 µL cass, 18": 320 µL; 1 µL cass, 30": 530 µL 10 µL cass, 18": 555 µL; 10 µL cass, 30": 920 µL
Dispense speed (high flow rate):	1 µL cass, 1 µL/well: 19s (96 wells) 55s (384 wells) 5 µL cass, 5 µL/well: 19s (96 wells), 58s (384 wells) 10 µL cass, 10 µL/well: 21s (96 wells), 66s (384 wells)

Physical Characteristics

Dimensions:	Base instrument: 17.19"W x 11.75" D x 8" H (43.51 x 29.21 x 20.32 cm)
Weight:	Base instrument: 19.5 lbs (8.8 Kg)
Power:	100 - 240 Volts AC. 50/60 Hz, 90 W max consumption
Connectivity:	Two USB ports: Protocol storage/transfer and for optional external mouse or keyboard

Regulatory

CE and TUV marked. RoHS compliant. IVD configurations are available.

Technical details are subject to change.



BioTek
BioTek Instruments, Inc.
Highland Park, P.O. Box 998
Winooski, Vermont 05404-0998, USA
Phone: 802-655-4040 • Toll-Free: 888-451-5171
Outside the USA: 802-655-4740
www.biotek.com