



The most sensitive
The most reliable



Helium Leak Detection

One Technology...Many applications

- CCI for container closure development
- CCI analysis at cryo storage conditions
- The most sensitive flow based leak test

Leak Detection Associates, a PTI company, is a global leading manufacturer of helium-based leak testing systems. For over 20 years, LDA has served the pharmaceutical, biotechnology, medical device and food packaging industries with solutions that encompass instrumentation, unique test fixtures, and all supporting services. LDA's offerings and industry experience have been used to establish global test standards (ASTM F2391-05) for these regulated industries, enabling package development, qualification, and routine testing to the highest standard. Preventative maintenance, service, and calibration programs are tailored to meet quality system standards and ensure proper instrument operation for years to come.


Why use Helium?

Helium leak detection (HeLD) is routinely used and widely accepted for applications that require the utmost leak sensitivity. Helium gas is an element with the second smallest atomic mass, is inert (non-reactive), is naturally rare in the atmosphere, and is non-flammable. Due to its low atomic mass this gas penetrates through the smallest cracks and micro-holes; being rare in the atmosphere eliminates background helium "noise" resulting in an extremely high signal-to-noise ratio thereby greatly increasing test sensitivity; being inert (non-reactive) ensures that the helium gas will not negatively affect the components being tested; and being non-flammable gives confidence that it is the safest tracer gas to use in a lab or production setting.

A custom-tuned spectrometer quantitatively measures the amount of helium escaping a package under vacuum. Leakage is quantified to levels well below Maximum Allowable Leakage Limits (MALL) of most all product-package systems, ideal for package development qualification and validation. Custom designed and easy to change test fixtures enable a wide range of applications. HeLD is currently used to fulfill test requirements and procedures in FDA guidance, EU Annex 1, ASTM F2391-05, and USP <1207>.

Benefits:

- Helium enables the discovery of extremely small microleaks that other leak testing methods cannot detect.
- Using a high vacuum technique, the leak test limits can be set down as low as 1×10^{-10} mbar-L/sec, a sensitivity level allowing unique comparisons between packaging components, materials, formats, and production parameters.



Cut-a-way view of the analyzer cell of a Helium Mass Spectrometer, showing how the incoming gas stream is separated. Only the Helium Ions follow a discrete path to the collection plate, where the signal is converted into a helium flow rate expressed in mBar-L/sec or atm-cc/sec.

SIMS 1915+

The Seal Integrity Monitoring System (SIMS) 1915+, is the optimal solution for helium-based leak detection for a variety of pharmaceutical and medical device product package systems. Common applications include vials, syringes, cartridges and blister cards. For parenteral products, common applications for helium include leak testing the rubber stopper on a vial or the plunger on a syringe or cartridge assembly.


Using helium as the tracer gas, packages can be quantitatively tested to levels far exceeding the vacuum bubble and dye penetration test methods. This quantitative approach allows direct comparison across various packaging materials and formats, production line settings and stability storage conditions, supporting the entire lifecycle.

This SIMS 1915+ will enable the quantitative analysis of packages at a sensitivity level as low as 1×10^{-10} mbar/L/sec. and provides relevant data sets in place of a simple pass/fail criteria while enabling testing to be performed at room temperature.

Each SIMS 1915+ Helium Leak Testing instrument is custom built to client specific standards and package configurations. We specialize in the engineering and development of custom test fixtures tailored to the component to be tested, which ensures precision and accuracy to meet your study goals, package configurations and quality monitoring needs.

SIMS 1915+ Leak Detection System includes the following standard components:

- Helium Leak Detector Module (HLDM): Oil Free, Production version in Console Frame Assembly, with Stainless Steel working surface.
- Vacuum Test Fixture Model (VTFM) - custom designed for use in testing package or container samples.
- Nitrogen Vent Supply system for the Helium Leak Detector Module (HLDM).
- Data Acquisition & Analysis Module running SIMS (ver. 2.0.5) 21 CFR Part 11 Compliant Leak Rate & Concentration Custom Software.
- External Calibrated Helium Leak Standard for System Calibration/Validation.
- Headspace Analyzer Module (HSAM) - Model VM-2 or UM-2A that includes three (3) Puncture Probes (two spares) for determination of Helium Concentration in the headspace of vial containers or blister card cavities. Includes Integral Calibration Chamber (and Calibrant Gas Flow Meter/Controller) provided for HSAM Calibration, and a Universal Holder to be used for various vial diameters.
- Dual Test Port Manifold (allows concurrent use of VTFM & HSAM on HLDM unit).
- Rotary Valve Box - for 3 Helium Calibration Gases (HSAM) & Nitrogen Gas (HLDM) Inlet Vent. Helium Filling Device specific for vial testing sample prep to enable vials to be filled with helium.



HELIUM FILL DEVICE

HEADSPACE ANALYZER MODULE (HSAM)



TEST CHAMBER FOR VIALS

SIMS 1915+ Options:

- LT 80 and LT 150 Low Temperature Add-on Systems for ultra-cold CCI testing using helium leak methodology.
- Various sizes of Vacuum Test Fixtures to accommodate all package types and sizes.
- Custom Test Fixtures to allow Helium leak testing for non-vacuum chamber applications.
- Additional External Helium Leak standards in various leak rate ranges.
- Helium Sniffer probe - for site specific determination of leak sites.
- Vacuum clamps, O-ring seals, chamber gaskets, and other spare parts.

SPECIFICATIONS

	SIMS 1915+
Test Method	Helium
Application	Testing and evaluation during package design and development, tooling qualification, production line setup and on-going product quality monitoring.
Package Type	Cold form blister cards, foil pouches, parenteral vials, cartridges, pre-filled syringes, bottles, combination product systems and medical device products.
Test Configuration	Console frame assembly with stainless steel working surface, including an articulating arm for mounting of computer monitor and keyboard.
Operator Interface	Customer supplied Window 10 (Pro) PC interface
Minimum Detectable Leak Rate	1×10^{-10} mbar L/sec.
CFR Security Capability	Data Acquisition & Analysis Module running SIMS (ver. 2.0.5) 21 CFR Part 11 Compliant Leak Rate & Concentration Custom Software
System Dimensions	25.5" (764mm) W x 30.25" (508mm) D x 38.5" (1153mm) H
Dimensions with Monitor arm	44.0" (111.8cm)/63.0" (160.0cm) Additional width and height
Weight	258 lbs. (117.3 kg.)
Power	100-240 VAC: 50/60 cycles
Options	Validation Qualification Package (IQ/OQ/PQ)
Test Method	ASTM F2391-05, Vacuum & Sniffer Modes Recommended in USP <1207> FDA Recognized Standard - Standard Test Method for Measuring Package and Seal Integrity Using Helium as the Tracer Gas
Backing Pump	Oil Free (No Exhaust) Backing Pump
Max Inlet Pressure	1.3 kPA
Operating Temperature	10-40° C

Low Temperature Add-on Test Systems

Complex drug and biologic formulations have resulted in life science companies continuous drive toward deeper cold storage in an effort to maintain product quality attributes. These products, often cell or gene therapies, or proteinaceous in nature, often require storage at temperatures below -20°C and are involved in storage and distribution environments incorporating dry ice (~-78.5°C), or even liquid nitrogen (~-200°C).

While the products demanding such intense cold storage may be complex, oftentimes, the package systems in which the products are placed are rather traditional in nature, such as a screw or crimp top vial. However, at these temperatures, many of the materials used in these package systems and responsible for maintaining package integrity are not typically assessed at these temperatures. When exposed to deep-cold or ultra-cold temperatures, physical changes to elastomeric components in particular can occur as materials reach or exceed their glass transition state, creating leaks at low temperatures that would otherwise not be observed while at room temperature. This type of leakage is typically observed at primary seal areas, such as that between an elastomeric closure and glass vial being used below -60C. Having a means to test container closure integrity while at these low temperatures enables manufacturers to gain insight into optimal package choice and design, as well as assembly parameters, to minimize leakage and demonstrate robust understanding of a package system's performance in accordance with USP <1207>.



LT 150 CHILLER FOR ULTRA COLD TEMP CCI TESTING ~-150° C.

SIMS 1915+ WITH LT 80 LOW-TEMPERATURE ADD-ON

To meet this market need for evaluating leaks at cold temperatures, LDA has created the LT80, Low Temperature -80°C, Add-On Test System for use with the LDA SIMS helium leak detectors. The LT80 system allows for concurrent temperature conditioning, temperature monitoring, and helium leak testing of packages approaching -80°C. The LT150 add-on low temperature unit enables sample testing close to -150° C. The LT Add-on units are available as upgrades to existing installed LDA SIMS units, or as a complete package.

APPLICATIONS




Vials

One of the most common package configurations for a variety of pharmaceutical and biotechnology products. Helium is the most ideal and sensitive method for component qualification for empty components as well as product filled vials.


Pre-Filled Syringes

The use of PFS systems has become more prevalent with the introduction of new and unique biological products. Helium test methods can be designed and developed that enable highly effective test programs for these unique package systems.



Bottles

Similar to vial package systems, bottles, whether composed of plastic materials or glass, continue to be one of the most widely used package systems and Helium remains the gold standard for quality testing purposes.



Blister Cards

Regardless of the cavity size and count, Helium remains the optimal test method for the qualification of the material components of CFF blister cards and has a long and trusted history for production line quality control verification.

Foil Pouches

The use of foil for pharmaceutical package systems remains in wide-spread use and Helium test methods continue to be a viable and highly sensitive approach to meeting strict regulatory requirements.



Combination Product Systems

Multi-chamber systems that require unique test requirements are ideally suited for the use of Helium leak testing to ensure all components meet the strict leak rate requirements.

Headquarters



8 Skyline Drive
Hawthorne, NY 10532

info@ptiusa.com
ptiusa.com

914.337.2005

Helium Center of Excellence



6638 Delilah Rd Suite F,
Egg Harbor Township, NJ 08234

info@ptiusa.com
heliumleak.com

609.415.2290
1.888.844.3759