



**AI-powered  
precision  
for colony  
counting**

Scan 3000 Ai



interscience  
Scan 5000 Ai



Scan Ai

Automatic AI colony counters  
& inhibition zone readers

interscience



## 20 years colony counting expertise

- Wide range of colony counters: manual, semi-auto, auto, real-time
- Dedicated AI R&D team since 2019
- More than 5000 **Scan**® equipment used worldwide every day
- Database of 1 million annotated images

SCIENTIFIC  
EQUIPMENT  
MANUFACTURER

**MADE IN FRANCE**



# A leap forward in counting performance

The **Scan® Ai** automatic colony counter automates and standardizes the counting of characteristic and non-characteristic colonies for microbiological analyses.


The power of Artificial Intelligence in the **Scan® Ai** automatic counter boosts Petri dish analysis with unrivalled speed and accuracy. **The resulting accuracy is 25% higher than that of a standard counter.**



Automatic colony counting in 1 second



Consistent accuracy up to 98%



400 plates counted per hour




Traceability of results

## Unrivalled precision with AI

Automatic counting ignores artifacts and counts on 100% of the plate. Artificial intelligence accurately detects and counts colonies, while distinguishing between different types of micro-organisms (bacteria, yeast, mold).

Counting improvements by AI based machine learning

Counting on 100% of the surface of the plate



Clear distinction between bacteria, yeasts and molds

Discrimination of artifacts (bubbles, marking, stickers...)

# The combination of powerful software and an ultra-precise counter

## Ultra HD optics

12.2 megapixel ultra HD camera<sup>(1)</sup>  
x 69 digital zoom

## Secure data

Works without an internet connection  
Data security and confidentiality  
Computer data stored locally



## Data integrity

Bi-directional connection  
Save and export results  
Double electronic signature  
Audit trail



IN ACCORDANCE WITH



(1) Only valid on Scan® 5000 Ai, please refer to technical specifications.



## AI colony counting

- 1000 colonies in 1 second
- Up to 98% counting accuracy
- AI learning from annotated images

## Inhibition zone reader

- Up to 16 antibiotics on a plate in 2 to 4 s<sup>(1)</sup>
- Measurement accuracy:  $\pm 0.1 \text{ mm}^{(1)}$
- Results categorized RIS



### Reflection- and shadow-free lighting

Diffusing White LED Dome lighting highlights colonies on the surface, pour plates and around the edges of the plate

### Wide reading range

Round Petri dishes from  $\varnothing 55$  to 150 mm<sup>(1)</sup>  
120 mm square Petri dishes<sup>(1)</sup>  
Compatible with multiple growth medias

### Sturdy

Robust frame in 304L stainless steel  
Impact-resistant glass  
3-year parts and labor warranty  
(after registration)

INTEGRATED DATABASE



# 25% greater accuracy than standard counters

Total count accuracy

73%

Classic counter

95%

Scan Ai

TBX positive colonies accuracy

73%

Classic counter

95%

Scan Ai

MRS accuracy

70%

Classic counter

95%

Scan Ai

VRBL accuracy

57%

Classic counter

89%


Scan Ai

## Applications


Artificial intelligence allows **Scan® Ai** to count a wide range of media, techniques and applications and different culture media for **QC**. The equipment can count Spiral® plates and Petrifilm™ with high precision and repeatability. Interpretation of chromogenic media is also available.

Integrated AI analysis of a wide range of different environments, techniques and growing media with accuracy and repeatability, including filtration membranes for liquid samples and **EM** contact plates. AI classification of different colony types, including bacteria, yeasts and molds, is also available.


**Scan® Ai** reads round plates with up to 12 antibiotic discs and square plates with up to 16 discs<sup>(1)</sup>




TBX




Air analysis  
on TSA




Discs  
on blood agar




TSA



Rhapsody




MRS




(1) Valid on the Scan® 5000 Ai model




# Comprehensive performance of pre-trained models





Square Petri dishes 120 mm<sup>(1)</sup>




Blood agar




Discs on agar



BAIRD PARKER



MILLIFLEX™ OASIS™



Mossel

The AI of **Scan® Ai** automatic counter has been trained on over a million dishes and annotated images from a wide range of industries.


We have integrated counting models on a large range of existing culture media:

- Total count (PCA/TSA)
- Coliforms
- Enterobacteriaceae (VRBL/VRBG)
- Escherichia Coli (TBX)
- Lactobacilles (MRS)
- Staphylococci (BAIRD-PARKER, BAIRD-PARKER RPF)
- Mossel
- Symphony
- Yeasts / molds
- GVPK
- Compass Bacillus Cereus
- Rhapsody
- SDA
- PDA
- Blood Agar
- MILLIFLEX™ OASIS™ (TSA, SDA, R2A).

**Scan® Ai** reads round Petri dishes up to 150 mm in diameter, and square 120 mm<sup>(1)</sup> plates.

It offers a wide choice of media and plates for greater flexibility:  
 Surface/ Pour / Mass / Spiral® / Circle  
 Petrifilm™, CompactDry™  
 chromogenic media,  
 MC-Media Pads™, Easy Plate™  
 Filtration membranes  
 MILLIFLEX™ OASIS™, Contact plates,  
 ATB on 120 mm square plates  
 ATB on blood agar

# Performance and safety with locked AI



## Continuous improvements in AI performance

Using artificial intelligence (AI) learning to count with convolutional neural networks (CNN) is transforming the field of colony counting.

When you choose **Scan® Ai**, you benefit from **high-performance, scalable automatic colony counting at the cutting edge of AI technology.**

## Enhanced safety with locked AI

The **Scan® Ai** system works autonomously, **without the need for an Internet connection.**

With locked Artificial Intelligence, **data is stored locally** for an even higher level of security. Local storage offers total control over information, keeping your data secure.



Update according to your technical  
and standards constraints



Keep control of your analyses  
Existing qualifications are retained




Evolve neural models with  
the analytical and matrix ecosystem





Software updates and AI models  
are available at your convenience



# Traceability and data integrity



**Scan® Ai** receives and returns information from the LIMS and/or dataLink™ and digitizes the count. **All your analysis data is traceable and auditable.**



**Results are validated by double electronic signature.** Thanks to account management, you have multi-level access to validate results.

This solution ensures data reliability and security in compliance with FDA guidelines, 21 CFR Part 11 and GMP (Good Manufacturing Practices) Annex 11.

## dataLink™, all your data on one plate



With dataLink™, traceability is guaranteed.

Set the label parameters which will be integrated on the sticker on the side of the plate.

# Scan<sup>®</sup> Ai range



**Scan<sup>®</sup> 3000 Ai**

Automatic AI colony counter & inhibition zone reader

Artificial intelligence

Ref. 434 300



**Scan<sup>®</sup> 5000 Ai**

Automatic AI colony counter & inhibition zone reader

Artificial Intelligence for large Petri dishes

Ref. 434 500

WEIGHTS AND DIMENSIONS	Scan <sup>®</sup> 3000 Ai	Scan <sup>®</sup> 5000 Ai
Dimensions (w x d x h)	32 x 32 x 44 cm	46.4 x 46.4 x 63.1 cm
Net weight	11.5 kg	25.15 kg
Box dimensions (w x d x h)	60 x 50 x 55 cm	80 x 60 x 63 cm
Gross weight	17.5 kg	32.50 kg

Supplied with 1 x 15 V power cord, 1 x USB cable, 3 validation plates, 1 x user manual, 1 x certificate of conformity / warranty card

Certified product

WEEE

2002/96/EC

RoHS

2002/95/EC






Products manufactured for INTERSCIENCE by Interlab, an ISO 9001-certified company

## Accessories



**Optical cleaner wipes**

Optical cleaner and 50 wipes

Ref. 435 702



**Barcode reader**

Barcode reader (1D / 2D)

Ref. 522 000



## IQOQ Services

In an environment where equipment precision and reliability are paramount, our application team is committed to providing you with **rigorous, customized qualification services**.

Our offers include a complete range of services, installation and initial commissioning to maintain the optimum performance of your equipment.

We understand the importance of each step in the process, and ensure that **each piece of equipment is qualified according to precise, validated protocols**.

# Technical specifications

	Scan® 3000 Ai	Scan® 5000 Ai
OVERVIEW		
Reference	434 300	434 500
Painted stainless steel shell	✓	✓
LIMS/SIL connection	✓	✓
USB connection	✓	✓
Available with dataLink™/dataLink™ pro traceability system	✓	✓
Counting on pour, surface, Spiral® and circle plated Petri dishes	✓	✓
Counting on chromogenic dishes	✓	✓
Counting on Petrifilm™, Compact Dry™, MC-Media Pads™, Easy Plate™, filtration membranes	✓	✓
Automatic counting	✓	✓
Inhibition zone reading	✓	✓
Counting on 100% of the Petri dish	✓	✓
AI-powered colony counting	✓	✓
Automatic detection of Petri dishes	-	✓
Counting on Petri dishes up to 150 mm	-	✓
COUNTING		
Counting	Automatic with manual control (adding/removing colonies)	
Automatic separation of clustered colonies	✓	✓
Creation of polygonal exclusion zones	✓	✓
Classification of bacteria, yeasts and molds	✓	✓
Counting time	Up to 1000 colonies per second	
Minimal size of colony	0.03 mm	
INHIBITION ZONE READING		
Antibiotic disc detection	Automatic with possibility to manually add or remove	
Automatic detection of antibiogram support	Disks (several brands simultaneously), wells, peni-cylinder (steel, plastic)	
Display resolution	± 0.1 mm	
Inhibition zone measurement accuracy	± 0.2 mm	± 0.1 mm
Number of antibiotic paper disks	Up to 7 antibiotics on a Ø 90 mm Petri dish	Up to 16 antibiotics on 120 mm square dish
Reading time	7 inhibition zone reading between 1 to 3 s	16 inhibition zone reading between 2 to 4 s max.
Interpretation system	CA-SFM Human health / EUCAST / CA-SFM Veterinary / CLSI (Clinical, Laboratory Standards Institute) / Customizable list	
SPECIFICATIONS		
Color camera	Ultra HD CMOS	
Lens	HD Japanese lens	
Zoom	x 69	
Resolution (megapixels)	5	12.2
White LED Lighting technology	White LED Dome indirect lighting	
LED Lighting system	Automatic with 7 combinations, top and/or bottom light, white or black background	
Petri dish size	Ø 55 mm - Ø 90 mm	Ø 55 mm to 150 mm round Petri dishes and 120 mm square Petri dishes
Color detection	4 colors on the same dish + 2 colors to exclude	
Languages	English, French, Japanese, Chinese, Russian, Spanish, German	
Voltage - Frequencies	100-240 V~ 50/60 Hz	
Warranty	3-year (after recording the warranty online)	
Spare parts availability	10 years	
In compliance with	21 CFR Part 11, ISO 7218 and AOAC 977.27	
TRACEABILITY		
USB Data export	Recountable session, Excel™, PDF report, jpg, png and bmp images	
Data security	Modified data traceability in conformity with 21 CFR part 11	
Results/traceability	Image / sample number / comments / date / time / antibiotic name / bacterial name / measured diameter / result categorized according to standards / minimum and maximum critical diameter	
PC MINIMUM REQUIREMENTS		
Operating sytem	Windows™ 10 or 11 (or higher)	
Processor	Intel i7, 2.8 GHz or higher (i9 or xeon gold)	
Graphic card	Nvidia RTX 3060 or 4050 or higher (Nvidia brand only)	
RAM	Minimum 16 GB required for use of Scan	
Equipments	Free USB 3.0 port	
Screen	1920 x 1080 pixels or higher	



# Plate & Count system®

The winning team  
from sample to result



**Plate & Count system®** enables **automatic dilution, plating and counting of colonies**.  
It's the perfect solution for the efficiency and traceability requirements of microbiological analysis.

- **75% savings in time and consumables guaranteed**
- **Traceability of results**
- **High-tech, made in France**

For more information on Plate & Count system, please refer to the brochure.

Your sales contact

## interscience

### PARIS

Phone: +33 (0)1 34 62 62 61 - Email: [info@interscience.com](mailto:info@interscience.com)

### FRANCFORT

Phone: +49 611 7238 7770 - Email: [sales.germany@interscience.com](mailto:sales.germany@interscience.com)

### BOSTON

Phone: +1 781 937 0007 - Email: [sales.usa@intersciencelab.com](mailto:sales.usa@intersciencelab.com)

### SHANGHAI

Phone: +86 (0)21-64739390 - Email: [sales.china@interscience.cn](mailto:sales.china@interscience.cn)

### SINGAPOUR

Phone: +65 6977 7232 - E-mail: [sales.asia@interscience.com](mailto:sales.asia@interscience.com)

### TOKYO

Phone: +81 3 6712 9715 - Email: [sales.japan@interscience.com](mailto:sales.japan@interscience.com)