

MIX-25/25P/30S Mini Mixer

Rotator, Mixer and Vortexer

Introduction

Mini Mixer is used for vibration, blend and stirring of sample tissue, cell, bacteria liquid and chemical reagent in life science and physical and chemical analysis field.

Feature

1. Mini and compact appearance, powerful function, sucker type under-chassis, super shockproof, ideal for high speed operation.
2. Multiple mixing modes, touch-operation, continuous operation and timing set.
3. Speed adjustable, and wide speed range, Stable low speed, powerful high speed.
4. Shaking blocks convenient for replacement, stable and reliable, sturdy and durable eccentric bearing designed.
5. Various tube holders are for optional (Suitable for Eppendorf Tubes).



Module fastener



MIX-30S



MIX-25/25P

Parameter

| Model | MIX-25 | MIX-25P | MIX-30S |
|-------------------------------|-----------------------------|-----------------------------|---------------------------------|
| Operation | knob + scale | knob + scale | key + digital Display |
| Shaking Orbit | 4mm | 3mm | 3mm |
| Shaking Mode | circle | circle | circle |
| Operation Mode | continuous or press running | continuous or press running | continuous running or short mix |
| Speed Range | 0~2500rpm | 0~2500rpm | 300~3000rpm |
| Time Range | / | / | 1s~999min |
| Voltage | AC 220V / AC 110V, 50/60Hz | AC 220V / AC 110V, 50/60Hz | DC24V |
| Power | 60W | 60W | 30W |
| Suitable Optional Accessories | MX-A and MX-B | all | all |
| Fuse | 250V, 1A, Φ 5x20 | 250V, 1A, Φ 5x20 | / |
| Dimension | W.132 x D.156 x H.150mm | W.132 x D.156 x H.150mm | W.132 x D.156 x H.150mm |
| Net Weight | 2.5kgs | 2.5kgs | 2.2kgs |