

# RH-24 3D Gyrotory Rocker

Shaker & Rocker

## Introduction

3D Gyrotory Rocker is an economical 3D gyrotory rocker with perfect 3-dimensional motion for gentle but thorough mixing. It is an ideal choice for mixing small sample. It suits for blood collection tubes and centrifuge tubes, applied for all blot related applications such as gel staining and destaining. The concave pad can firmly keep the tubes.

## Feature

1. Gentle but thorough mixing.
2. Suits for blood collection tubes and centrifuge tubes.
3. Fixed speed and tilt angle.
4. Safe for use in cold rooms and incubators.

## Parameter

Model	RH-24
Speed Range	20rpm (fixed)
Shaking Orbit	3D / 20 (fixed)
Max. Load	0.8kgs
Voltage	AC 220V/AC 110V, 50/60Hz
Power	45W
Dimension	W.240 x D.170 x H.150mm
Net Weight	1.6kgs



## Remarks

RH-24 is available for 1.5mlx60, 3mlx30, 7mlx21, 15ml x17.  
Also suitable for use in cold closet.

